


# Apera LabSen® Professional pH Electrodes

Apera Instruments LabSen® Professional pH Electrodes are backed by proprietary sensor technologies and components from Switzerland, designed for a wide variety of applications in scientific research and product quality control with more than 30 models.



## Glass Membrane

- Glass membrane is the most important part of pH electrodes.
- LabSen® pH electrodes are equipped with 4 types of glass membrane to meet needs in various applications: S membrane, H membrane, HF membrane, and PHY membrane.
  - LabSen® glass membrane is highly resistant to general impact (totally differentiates from the traditional fragile glass bulb membrane).
  - LabSen® glass membrane with different shapes are shown as below:



Hemispherical Cylindrical Slim Spear Flat Conical

## Junction

- Junction is the electrolyte interface between reference system and the test samples. LabSen® electrode adopts the following junctions:
- Diaphragm – the most frequently used junction type, easy to get blocked by protein-containing or suspension solutions.
  - Pore without Diaphragm – it is used with solid electrolyte, no clogging, and maintenance-free.
  - Moveable sleeve – easy to clean, suitable for suspension, emulsion, low ion-concentration solution and non-aqueous solution. The infiltration rate of electrolyte is determined by the tightness of the sleeve during installation.
  - PTFE – a type of Teflon material with multi pores, hard to be contaminated.

## Inner Solution

The inner solution of LabSen® electrode is in a distinctive dark blue color. With a special gel treatment, inner solution does not flow and will not cause bubble. The electrode can even work upside down.

## Reference System

Besides the routine Ag/AgCl reference electrode, LabSen® Electrodes are more likely to adopt the Long-Life reference system and Silver Ion Trap reference electrode.



Long-Life Reference System

Long-Life reference system is composed of a glass tube, AgCl, and reference silver wire. The top end of the slim glass tube is stuffed with cotton, which will prevent reaction between AgCl and electrolyte when temperature changes. It improves the stability of reference electrode and service life.

Material	Lead-Free Glass
Temp. Probe	No
Junction	Ceramic
Reference	Long Life
Electrolyte	3M KCl
Length	120mm
Diameter	12mm
Membrane Shape	Half Ball

### LabSen® 211 Routine pH Electrode



Range: 0 to 14 pH  
Temp.: -5 to 100 °C  
(23 to 212°F)  
Connector: BNC/1m cable

#### Application

- Suitable for routine lab tests in regular aqueous solutions

#### Features

- Shock-proof glass membrane
- Long-Life reference

Material	Lead-Free Glass
Temp. Probe	No
Junction	Movable Sleeve
Reference	Long Life
Electrolyte	3M KCl
Length	130mm
Diameter	12mm
Membrane Shape	Cylindrical

### LabSen® 221 Precise pH Electrode



Range: 0 to 14 pH  
Temp.: -5 to 80 °C  
(23 to 176°F)  
Connector: BNC/1m cable

#### Application

- Suitable for low ion-concentration liquid and viscous samples.

#### Features

- Movable sleeve helps generate quick and stable readings in low ion-concentration samples.

Material	Lead-Free Glass
Temp. Probe	Probe
Junction	Ceramic
Reference	Long Life
Electrolyte	3M KCl
Length	150mm
Diameter	12-6mm
Membrane Shape	Cylindrical

### LabSen® 241-6 Semi-micro pH Electrode



Range: 0~14 pH  
Temp.: 0 to 100 °C  
(32 to 212°F)  
Connector: BNC/1m cable

#### Application


- Suitable for small volume liquid ( $\geq 0.2\text{ml}$ ) and test tube samples.

#### Features

- 6mm measuring tip
- Long-Life reference system

Material	Lead-Free Glass
Temp. Probe	No
Junction	Ceramic
Reference	Long Life
Electrolyte	3M KCl
Length	130mm
Diameter	12-3mm
Membrane Shape	Cylindrical

### LabSen® 241-3 Micro pH Electrode



Range: 0 to 14 pH  
Temp.: 0 to 100 °C  
(32 to 212°F)  
Connector: BNC/1m cable

#### Application

- Suitable for micro samples ( $\geq 30\mu\text{L}$ ), or testing in centrifuge tubes and NMR tubes.

#### Features

- 3mm measuring tip is  $\Phi 3 \times 70\text{mm}$ .
- Long-Life reference system

Material	Lead-Free Glass
Temp. Probe	No
Junction	Ceramic +Single Pore
Reference	Long Life
Electrolyte	Polymer
Length	100mm
Diameter	12-6mm
Membrane Shape	Spear

### LabSen® 251 Glass Spear pH Electrode



Range: 0 to 14 pH  
Temp.: 0 to 80 °C  
(32 to 176°F)  
Connector: BNC/1m cable

#### Application

- Suitable for soft solid and semi-solid food, e.g. cheese, fruit, vegetables, sushi rice, etc.

#### Features

- Spear tip for direct soil pH measurements
- Solid Polymer electrolyte, maintenance-free

Material	PVC
Temp. Probe	No
Junction	Ceramic +Single Pore
Reference	Long Life
Electrolyte	Polymer
Length	90mm
Diameter	15-6mm
Membrane Shape	Spears

### LabSen® 551 Spear pH Electrode



Range: 0 to 14 pH  
Temp.: 0 to 60°C  
(32 to 140°F)  
Connector: BNC/1m Cable

- Application**
- Suitable for solid and semi-solid medium e.g. soil
- Features**
- Spear tip for direct soil pH measurements
  - Solid Polymer electrolyte, maintenance-free

Material	POM
Temp. Probe	No
Junction	PTFE
Reference	Long Life
Electrolyte	3M KCl
Length	110mm
Diameter	12mm
Membrane Shape	Flat

### LabSen® 371 Flat pH Electrode



Range: 0 to 14 pH  
Temp.: 0 to 80°C  
(32 to 176°F)  
Connector: BNC/1m Cable

- Application**
- Suitable for flat surface measurement such as fabrics, skin, paper, carpet, etc.
- Features**
- Flat glass membrane, fast response

Material	Lead-free Glass
Temp. Probe	No
Junction	Movable Sleeve
Reference	Silver Ion Trap
Electrolyte	3M KCl
Length	130mm
Diameter	12mm
Membrane Shape	Cylindrical

### LabSen® 801 Purified Water pH Electrode



Range: 1 to 11 pH  
Temp.: 0 to 80°C  
(32 to 176°F)  
Connector: BNC/1m Cable

- Application**
- Suitable for purified water measurement e.g. RO water, distilled water, deionized water etc.
- Features**
- Movable sleeve helps generate quick and stable measurements in low ion-concentration liquid.
  - Silver ion trap reference system

Material	Lead-Free Glass
Temp. Probe	No
Junction	Ceramic
Reference	Silver Ion Trap
Electrolyte	3M KCl
Length	120mm
Diameter	12mm
Membrane Shape	Half Ball

### LabSen® 831 HF pH Electrode



Range: 0 to 11 pH  
Temp.: 0 to 100°C  
(32 to 212°F)  
Connector: BNC/1m Cable

- Application**
- Suitable for strong acid solutions and solutions containing Hydrofluoric acid ( $\geq 3\text{pH}$ ).
- Features**
- Special HF glass membrane, HF corrosion resistance.
  - Silver ion trap reference system, preventing the junction from being contaminated by silver sulfide or protein.

Material	Lead-Free Glass
Temp. Probe	No
Junction	Ceramic *3
Reference	Long Life
Electrolyte	Protelyte
Length	120mm
Diameter	12-6mm
Membrane Shape	Conical

### LabSen® 851-1 Viscous pH Electrode



Range: 0 to 14 pH  
Temp.: 0 to 100°C  
(32 to 212°F)  
Connector: BNC/1m Cable

- Application**
- Suitable for cosmetics and viscous samples.
- Features**
- 3 ceramic diaphragm
  - Protelyte reference solution
  - Silver ion trap reference system